



# MEDICAL FACE MASK

## EARLOOPS • ANTI-GLARE / ANTI-FOG VISOR

### ANTI-FOG FOAM • SOFT • HYPOALLERGENIC

### WHITE • ASTM LEVEL 3



**312251EMV**

Tested against the ASTM F2100 Standards, this ASTM Level 3 mask offers high filtration efficiency and is designed to withstand intense spray impact in areas such as trauma, ICU, and the ER. Its attached wrap-around visor, with built in anti-fog properties, and an anti-fog foam strip provide clinicians with optical clarity and enhanced protective features. The premium inner layer — a soft, irritant-free material — makes this the perfect mask for individuals with skin sensitivities.

BFE	PFE	DELTA P $\Delta P(\text{mmHg}/\text{cm}^2)$	FLUID RESISTANCE	FLAME SPREAD
> 98%	> 85%	< 6.0	160 mm Hg	Class 1
MASKS/BOX		BOXES/CASE		MASKS/CASE
25		4		100

**① Soft Irritant Free Inner Layer**

- Inner layer has a silky finish that is delicate yet comfortable, making it ideal for those with skin sensitivities.

**② High Quality Nosepiece with Ultra Soft Foam**

- Aluminum nosepiece and soft foam strip work together to form smoothly across the nose and cheekbones.
- Remains securely in place to create a custom and secure seal providing enhanced resistance to fogging of eyewear and a cushion-like feel to the mask.

**③ Anti-Fog Visor with Wrap-Around Design**

- Visor design provides clinicians with a clear line of sight paired while protecting from direct and peripheral protection from spray of fluids.

