



MEDICAL FACE MASK

EARLOOPS • ANTI-FOG VISOR

ANTI-FOG FOAM • OPAL • ASTM LEVEL 3



313141EMV

Tested against the ASTM F2100 Standards, this ASTM Level 3 mask offers high filtration efficiency and is designed to withstand intense spray impact in areas such as trauma, ICU, and the ER. Its attached wrap-around visor with built in anti-fog properties and an anti-fog foam strip provides clinicians with optical clarity and enhanced protective features.

BFE	PFE	DELTA P $\Delta P(\text{mmH}_2\text{O}/\text{cm}^2)$	FLUID RESISTANCE	FLAME SPREAD
> 98%	> 85%	< 6.0	160 mm Hg	Class 1
MASKS/BOX		BOXES/CASE		MASKS/CASE
25		4		100

- 1

**Ultra-soft Earloop**
  - Earloop material is soft and smooth for maximum comfort.
  - Round design makes this mask one of the most comfortable options available.
- 2

**High Quality Nosepiece with Ultra Soft Foam**
  - Aluminum nosepiece and soft foam strip work together to form smoothly across the nose and cheekbones.
  - Remains securely in place to create a custom and secure seal providing enhanced resistance to fogging of eyewear and a cushion-like feel to the mask.
- 3

**Anti-Fog Visor with Wrap-Around Design**
  - Visor design provides clinicians with a clear line of sight paired while protecting from direct and peripheral protection from spray of fluids.

