



MEDICAL FACE MASK

EARLOOPS • ANTI-FOG FOAM • OPAL

ASTM LEVEL 3



31314EM

This ASTM Level 3 mask offers the highest levels of fluid resistance and filtration efficiency with the added benefit of an anti-fog foam strip for enhanced visibility when worn with eyewear. Clinicians reach for this mask when performing medical procedures with a high risk of fluid splash in areas including the ER, ICU, and trauma.

BFE	PFE	DELTA P $\Delta P(\text{mmH}_2\text{O}/\text{cm}^2)$	FLUID RESISTANCE	FLAME SPREAD
> 98%	> 85%	< 6.0	160 mm Hg	Class 1

MASKS/BOX	BOXES/CASE	MASKS/CASE
50	10	500

- ① Ultra-soft Earloop**
 - Earloop material is soft and smooth for maximum comfort.
 - Round design makes this mask one of the most comfortable options available.
- ② High Quality Nosepiece with Ultra Soft Foam**
 - Aluminum nosepiece and soft foam strip work together to form smoothly across the nose and cheekbones.
 - Remains securely in place to create a custom and secure seal providing enhanced resistance to fogging of eyewear and a cushion-like feel to the mask.
- ③ Smooth Naturally Breathable Inner Layer**
 - Provides a light and cool feel to the mask, reducing chances of irritation throughout extended wear.
 - Inner layer maintains its structure, allowing the mouth to feel open and airy.

