



MEDICAL FACE MASK

TIES • ANTI-FOG VISOR • ANTI-FOG FOAM

OPAL • ASTM LEVEL 3



313241TMV

Tested against the ASTM F2100 Standards, this ASTM Level 3 mask provides high filtration efficiency and is designed to withstand intense spray impact throughout surgical procedures in areas such as trauma and the ER. Its attached wrap-around visor with built in anti-fog properties and an anti-fog foam strip provides clinicians with optical clarity and enhanced protective features.

BFE	PFE	DELTA P ΔP(mmH ₂ O/cm ²)	FLUID RESISTANCE	FLAME SPREAD
> 98%	> 85%	< 6.0	160 mm Hg	Class 1
MASKS/BOX		BOXES/CASE		MASKS/CASE
25		4		100

① **Premium Tie Material and Design**

- Tie material that is comfortable on the skin throughout extended wear.
- Long ties create a flexible fit while tying over top of a bouffant or surgical cap.

② **High Quality Nosepiece with Ultra Soft Foam**

- Aluminum nosepiece and soft foam strip work together to form smoothly across the nose and cheekbones.
- Remains securely in place to create a custom and secure seal providing enhanced resistance to fogging of eyewear and a cushion-like feel to the mask.

③ **Anti-Fog Visor with Wrap-Around Design**

- Visor design provides clinicians with a clear line of sight paired while protecting from direct and peripheral protection from spray of fluids.

