



MEDICAL FACE MASK

TIES • ANTI-GLARE / ANTI-FOG VISOR

ANTI-FOG FOAM • OPAL • ASTM LEVEL 3



313242TMV

Created using only premium materials, this ASTM Level 3 mask passes the highest standards for fluid resistance and filtration efficiency. Its attached wrap-around visor with anti-fog and anti-glare properties and an anti-fog foam strip ensures a clear line of sight, making this a go-to mask for even the most laborious surgical procedures.

BFE	PFE	DELTA P $\Delta P(\text{mmH}_2\text{O}/\text{cm}^2)$	FLUID RESISTANCE	FLAME SPREAD
> 98%	> 85%	< 6.0	160 mm Hg	Class 1
MASKS/BOX		BOXES/CASE		MASKS/CASE
25		4		100

① Premium Tie Material and Design

- Tie material that is comfortable on the skin throughout extended wear.
- Long ties create a flexible fit while tying over top of a bouffant or surgical cap.

② High Quality Nosepiece with Ultra Soft Foam

- Aluminum nosepiece and soft foam strip work together to form smoothly across the nose and cheekbones.
- Remains securely in place to create a custom and secure seal providing enhanced resistance to fogging of eyewear and a cushion-like feel to the mask.

③ Premium Anti-Glare / Anti-Fog Visor

- Built in anti-glare and anti-fog treatments provide better clarity, unencumbered by OR and other hospital lights.
- Wrap-around design provides clinicians with a clear line of sight while protecting from direct and peripheral spray of fluids.



Not made with
natural rubber latex



VIVIT Supply LLC

200 Continental Dr., Suite 401
Newark, DE 19713 USA

W vivitsupply.com

V004