



# OVERHEAD SPUNDBOND GOWN THUMBLOOPS WITH ELASTIC CUFF YELLOW



## 130Y UNIVERSAL / 131Y XL

Engineered to provide ultimate comfort, this gown is a go-to when fluid contact is not expected. Made from strong, comfortable and breathable polypropylene material combined with the open back design allows the wearer to stay cool and comfortable on the job.

Impact Penetration	Hydrostatic Pressure	Weight	Fluid Resistance	Breathability
-	-	Light		

SKU#	SIZE	GOWNS/BAG	BAGS/CASE	GOWNS/CASE
130Y	UNIVERSAL	10	10	100
131Y	EXTRA LARGE	10	10	100

### ① “Clean Release” Overhead Neck Closure

- Overhead neck closure promotes compliance and increases donning & doffing speed
- Removal of a shoulder seam creates a proper release at the back neck closure removing risk of contamination

### ② Soft & Breathable Spunbond

- Extra soft lightweight single layer of polypropylene for increased comfort
- Maximized breathability throughout the gown

### ③ Aligned Thumbloop with Welded Elastic Cuff

- Properly aligned thumbloops increases donning speed
- Patented welded elastic technology eliminates pinholes, reducing risk of contamination

