

OVERHEAD SMS GOWN

THUMBLOOPS WITH ELASTIC CUFF

YELLOW



230Y UNIVERSAL / 231Y XL

Engineered to provide ultimate comfort, this gown is a go-to when fluid contact is not expected. Strong breathable SMS material combined with the open back design allows the wearer to stay cool and comfortable on the job.

Impact Penetration	Hydrostatic Pressure	Weight	Fluid Resistance	Breathability
-	-	Light		

SKU#	SIZE	GOWNS/BAG	BAGS/CASE	GOWNS/CASE
230Y	UNIVERSAL	10	10	100
231Y	EXTRA LARGE	10	10	100

- ① **“Clean Release” Overhead Neck Closure**
 - Overhead neck closure promotes compliance and increases donning & doffing speed
 - Removal of a shoulder seam creates a proper release at the back neck closure removing risk of contamination
- ② **Soft & Breathable SMS**
 - Specially designed material maximizes breathability and comfort
- ③ **Aligned Thumbloop with Welded Elastic Cuff**
 - Properly aligned thumbloops increases donning speed
 - Patented welded elastic technology eliminates pinholes, reducing risk of contamination

