



OVERHEAD AAMI LEVEL 2 ISOLATION GOWN • OPEN BACK THUMBLOOP WITH ELASTIC CUFF • YELLOW



240Y UNIVERSAL / 241Y XL

Breathable SMS material utilized in this gown promotes strength and fluid resistance, ideally suited for tasks involving minimal to low risk of fluid exposure. This gown comes equipped with a convenient overhead neck closure and thumbloops with elastic cuffs to ensure reliable protection from exposure.

Impact Penetration	Hydrostatic Pressure	Weight	Fluid Resistance	Breathability
≤1.0 g	≥20 cm	Medium		

SKU#	SIZE	GOWNS/BAG	BAGS/CASE	GOWNS/CASE
240Y	UNIVERSAL	10	10	100
241Y	EXTRA LARGE	10	10	100

① “Clean Release” Overhead Neck Closure

- Overhead neck closure promotes compliance and increases donning & doffing speed
- Removal of a shoulder seam creates a proper release at the back neck closure removing risk of contamination

② Soft & Breathable SMS

- Specially designed material allows for the perfect mixture of fluid protection and comfort
- Maximizes breathability throughout the gown

③ Aligned Thumbloop with Welded Elastic Cuff

- Properly aligned thumbloops increases donning speed
- Patented welded elastic technology eliminates pinholes, reducing risk of contamination



Not made with
natural rubber latex

