

# POLY-COATED GOWN

## NECK TIES • KNIT CUFF



### YELLOW



#### 320Y UNIVERSAL / 321Y XL

Engineered to provide the ultimate balance between protection and comfort, this gown is a go-to in decontamination departments. Strong poly-coated spunbond material achieves protection from penetration of potentially harmful fluids and the neck ties allow adjustability for the wearer to stay cool and comfortable on the job.

Raw Material Performance meets Impact Penetration and Hydrostatic Pressure requirements at AAMI PB70 Level 3:

Impact Penetration	Hydrostatic Pressure	Weight	Fluid Resistance	Breathability
≤1.0 g	≥50 cm	Heavy		

Made with highly fluid resistant materials. Gown is not PB70 rated.

SKU#	SIZE	GOWNS/BAG	BAGS/CASE	GOWNS/CASE
320Y	UNIVERSAL	10	10	100
321Y	EXTRA LARGE	10	10	100

#### ① Neck Tie Closure

- Fastens securely at both the neck and waist for a snug fit
- Easy to adjust fit for maximum flexibility ensuring the highest level of coverage

#### ② Strong & Durable Poly-Coated Spunbond

- Outer layer of polyethylene film provides protection against fluids
- Inner layer of polypropylene allows for increased comfort

#### ③ Extra Soft Knit Cuff

- Material provides maximum comfort during periods of extended use.
- Snug fit to keep cuff in place underneath glove allowing for increased protection

