



OVERHEAD POLY-COATED GOWN THUMBLOOPS WITH ELASTIC CUFF YELLOW



330Y UNIVERSAL / 331Y XL

Engineered to provide the ultimate balance between protection and comfort, this gown is a go-to in decontamination departments. Strong poly-coated spunbond material provides protection from fluids and the open back design allows the wearer to stay cool and comfortable on the job.

Raw Material Performance meets Impact Penetration and Hydrostatic Pressure requirements at AAMI PB70 Level 3:

Impact Penetration	Hydrostatic Pressure	Weight	Fluid Resistance	Breathability
≤1.0 g	≥50 cm	Heavy	🔴🔴🔴	🌊🌊🌊🌊

Made with highly fluid resistant materials. Gown is not PB70 rated.

SKU#	SIZE	GOWNS/BAG	BAGS/CASE	GOWNS/CASE
330Y	UNIVERSAL	10	10	100
331Y	EXTRA LARGE	10	10	100

① “Clean Release” Overhead Neck Closure

- Overhead neck closure promotes compliance and increases donning & doffing speed
- Removal of a shoulder seam creates a proper release at the back neck closure removing risk of contamination

② Strong & Durable Poly-Coated Spunbond

- Outer layer of polyethylene film provides highest protection against fluids
- Inner layer of polypropylene allows for increased comfort

③ Aligned Thumbloop with Welded Elastic Cuff

- Properly aligned thumbloops increases donning speed
- Patented welded elastic technology eliminates pinholes, reducing risk of contamination

