



# OVERHEAD AAMI LEVEL 4 GOWN OPEN BACK • KNIT CUFF • BLUE



## 436B UNIVERSAL / 437B XL

Designed to deliver the perfect balance of protection and comfort, this gown is a top choice for decontamination departments. Made from strong, fluid resistant poly-coated spunbond material, this open back gown comes equipped with a convenient overhead neck closure and knit cuffs to ensure maximum comfort.

ASTM F1671	Weight	Fluid Resistance	Breathability
Pass	Heavy		

Variation in lot results exist. All results meet or exceed AAMI PB70 requirements.

SKU#	SIZE	GOWNS/BAG	BAGS/CASE	GOWNS/CASE
436B	UNIVERSAL	10	10	100
437B	EXTRA LARGE	10	10	100

### ① “Clean Release” Overhead Neck Closure

- Overhead neck closure promotes compliance and increases donning & doffing speed
- Removal of a shoulder seam creates a proper release at the back neck closure reducing the risk of contamination

### ② AAMI 4 Rated Material and Seams

- AAMI 4 Level 4-rated material and seams provide the highest level of AAMI rated protection
- Outer layer of polyethylene film provides highest protection against fluids
- Inner layer of polypropylene allows for increased comfort

### ③ Extra Soft Knit Cuff

- Material provides maximum comfort during periods of extended use
- Snug fit to keep cuff in place underneath glove allowing for increased protection



Not made with  
natural rubber latex



VIVIT Supply LLC  
200 Continental Dr., Suite 401  
Newark, DE 19713 USA

W vivitsupply.com

V001