



OVERHEAD AAMI LEVEL 2 ISOLATION GOWN • THUMBLOOP WITH ELASTIC CUFF • YELLOW



820Y UNIVERSAL / 821Y XL

Breathable SMS material utilized in this gown promotes strength and fluid resistance, ideally suited for tasks involving minimal to low risk of fluid exposure. This gown comes equipped with a convenient overhead neck closure and thumbloops with elastic cuffs to ensure reliable protection from exposure.

Impact Penetration	Hydrostatic Pressure	Weight	Fluid Resistance	Breathability
≤1.0 g	≥20 cm	Medium		

Variation in lot results exist. All results meet or exceed AAMI PB70 requirements.

SKU#	SIZE	GOWNS/BAG	BAGS/CASE	GOWNS/CASE
820Y	UNIVERSAL	10	10	100
821Y	EXTRA LARGE	10	10	100

- 1

“Clean Release” Overhead Neck Closure
 - Overhead neck closure promotes compliance and increases donning & doffing speed
 - Removal of a shoulder seam creates a proper release at the back neck closure removing risk of contamination
- 2

Soft & Breathable SMS
 - Specially designed material allows for the perfect mixture of fluid protection and comfort
 - Maximizes breathability throughout the gown
- 3

Aligned Thumbloop with Welded Elastic Cuff
 - Properly aligned thumbloops increases donning speed
 - Patented welded elastic technology eliminates pinholes, reducing risk of contamination

