



OVERHEAD AAMI LEVEL 3 ISOLATION GOWN • THUMBLOOP WITH ELASTIC CUFF • YELLOW



830Y UNIVERSAL / 831Y XL

Breathable SMS material utilized in this gown promotes strength and fluid resistance, ideally suited for tasks involving moderate risk of fluid exposure. This gown comes equipped with a convenient overhead neck closure and thumbloops with elastic cuffs to ensure reliable protection from exposure.

Impact Penetration	Hydrostatic Pressure	Weight	Fluid Resistance	Breathability
≤1.0 g	≥50 cm	Heavy		

Variation in lot results exist. All results meet or exceed AAMI PB70 requirements.

SKU#	SIZE	GOWNS/BAG	BAGS/CASE	GOWNS/CASE
830Y	UNIVERSAL	10	10	100
831Y	EXTRA LARGE	10	10	100

① “Clean Release” Overhead Neck Closure

- Overhead neck closure promotes compliance and increases donning & doffing speed
- Removal of a shoulder seam creates a proper release at the back neck closure removing risk of contamination

② Soft & Breathable SMS

- Specially designed material allows for the perfect mixture of fluid protection and comfort
- Maximizes breathability throughout the gown

③ Aligned Thumbloop with Welded Elastic Cuff

- Properly aligned thumbloops increases donning speed
- Patented welded elastic technology eliminates pinholes, reducing risk of contamination



Not made with
natural rubber latex

