



AAMI LEVEL 2 ISOLATION GOWN

NECK TIES • THUMBLOOP WITH ELASTIC CUFF • YELLOW



840Y UNIVERSAL / 841Y XL

Breathable SMS material utilized in this gown promotes strength and fluid resistance, ideally suited for tasks involving minimal to low risk of fluid exposure. This gown comes equipped with neck ties for adjustability and thumbloops with elastic cuffs to enhance protection from exposure.

Impact Penetration	Hydrostatic Pressure	Weight	Fluid Resistance	Breathability
≤1.0 g	≥20 cm	Medium		

Variation in lot results exist. All results meet or exceed AAMI PB70 requirements.

SKU#	SIZE	GOWNS/BAG	BAGS/CASE	GOWNS/CASE
840Y	UNIVERSAL	10	10	100
841Y	EXTRA LARGE	10	10	100

① Neck Tie Closure

- Fastens securely at both the neck and waist for a snug fit
- Easy to adjust fit for maximum flexibility ensuring the highest level of coverage

② Soft & Breathable SMS

- Specially designed material allows for the perfect mixture of fluid protection and comfort
- Maximizes breathability throughout the gown

③ Aligned Thumbloop with Welded Elastic Cuff

- Properly aligned thumbloops increases donning speed
- Patented welded elastic technology eliminates pinholes, reducing risk of contamination

