



AAMI LEVEL 3 ISOLATION GOWN

NECK TIES • THUMBLOOP WITH ELASTIC CUFF • YELLOW



850Y UNIVERSAL / 851Y XL

Breathable SMS material utilized in this gown promotes strength and fluid resistance, ideally suited for tasks involving moderate risk of fluid exposure. This gown comes equipped with neck ties for adjustability and thumbloops with elastic cuffs to enhance protection from exposure.

Impact Penetration	Hydrostatic Pressure	Weight	Fluid Resistance	Breathability
≤1.0 g	≥50 cm	Heavy		

Variation in lot results exist. All results meet or exceed AAMI PB70 requirements.

SKU#	SIZE	GOWNS/BAG	BAGS/CASE	GOWNS/CASE
850Y	UNIVERSAL	10	10	100
851Y	EXTRA LARGE	10	10	100

① Neck Tie Closure

- Fastens securely at both the neck and waist for a snug fit
- Easy to adjust fit for maximum flexibility ensuring the highest level of coverage

② Soft & Breathable SMS

- Specially designed material allows for the perfect mixture of fluid protection and comfort
- Maximizes breathability throughout the gown

③ Aligned Thumbloop with Welded Elastic Cuff

- Properly aligned thumbloops increases donning speed
- Patented welded elastic technology eliminates pinholes, reducing risk of contamination



Not made with natural rubber latex

