



AAMI LEVEL 2 ISOLATION GOWN • TAPE TAB ELASTIC CUFF • YELLOW



860Y UNIVERSAL / 861Y XL

Breathable SMS material utilized in this gown promotes strength and fluid resistance, ideally suited for tasks involving minimal to low risk of fluid exposure. This gown comes equipped with a convenient tape tab neck closure and elastic cuffs to ensure reliable protection from exposure.

Impact Penetration	Hydrostatic Pressure	Weight	Fluid Resistance	Breathability
≤1.0 g	≥20 cm	Medium		

Variation in lot results exist. All results meet or exceed AAMI PB70 requirements.

SKU#	SIZE	GOWNS/BAG	BAGS/CASE	GOWNS/CASE
860Y	UNIVERSAL	10	10	100
861Y	EXTRA LARGE	10	10	100

① Tape Tab Neck Closure

- Easy to adjust fit for maximum flexibility ensuring the highest level of coverage

② Soft & Breathable SMS

- Specially designed material allows for the perfect mixture of fluid protection and comfort
- Maximizes breathability throughout the gown

③ Welded Elastic Cuff

- Patented welded elastic technology eliminates pinholes, reducing risk of contamination

