

# AAMI LEVEL 3 ISOLATION GOWN • TAPE TAB ELASTIC CUFF • YELLOW



## 870Y UNIVERSAL / 871Y XL

Breathable SMS material utilized in this gown promotes strength and fluid resistance, ideally suited for tasks involving moderate risk of fluid exposure. This gown comes equipped with a convenient tape tab neck closure and elastic cuffs to ensure reliable protection from exposure.

Impact Penetration	Hydrostatic Pressure	Weight	Fluid Resistance	Breathability
≤1.0 g	≥50 cm	Heavy		

Variation in lot results exist. All results meet or exceed AAMI PB70 requirements.

SKU#	SIZE	GOWNS/BAG	BAGS/CASE	GOWNS/CASE
870Y	UNIVERSAL	10	10	100
871Y	EXTRA LARGE	10	10	100

### ① Tape Tab Neck Closure

- Easy to adjust fit for maximum flexibility ensuring the highest level of coverage

### ② Soft & Breathable SMS

- Specially designed material allows for the perfect mixture of fluid protection and comfort
- Maximizes breathability throughout the gown

### ③ Welded Elastic Cuff

- Properly aligned thumbloops increases donning speed
- Patented welded elastic technology eliminates pinholes, reducing risk of contamination

