



CHEMOTHERAPY DRUG TESTED AAMI LEVEL 2 ISOLATION GOWN KNIT CUFF • PURPLE



920P UNIVERSAL / 921P XL

VIVIT Supply Chemotherapy Drug Tested Isolation Gowns are designed to comply with all USP <800> requirements. Gown seams have been tested using ASTM F739 to give clinical staff the protection they need when handling or administering harsh chemotherapy drugs. This gown is ideally suited for tasks involving minimal to low risk of fluid exposure.

Impact Penetration	Hydrostatic Pressure	Weight	Fluid Resistance	Breathability
≤1.0 g	≥20 cm	Heavy		

Variation in lot results exist. All results meet or exceed AAMI PB70 requirements.

SKU#	SIZE	GOWNS/BAG	BAGS/CASE	GOWNS/CASE
920P	UNIVERSAL	10	5	50
921P	EXTRA LARGE	10	5	50

① Neck Closure — HOOK & LOOP

- Easy to adjust fit for maximum flexibility ensuring the highest level of coverage

② Material — Specialized SMS/PE

- High quality SMS/PE material, ultrasonically welded, ensuring maximum protection against harsh chemotherapy drugs
- Unique purple color to easily identify gown

③ Extra Soft Knit Cuff

- Material provides maximum comfort during periods of extended use.
- Snug fit to keep cuff in place underneath glove allowing for increased protection

