

OVERHEAD FILM GOWN

THUMBLOOPS • BLUE



970B UNIVERSAL / 971B XL

Made of a single, lightweight layer of protective polyethylene film which provides a strong level of protection from fluids. The open back design provides the breathability to keep the wearer cool and comfortable on the job.

Raw Material Performance meets Impact Penetration and Hydrostatic Pressure requirements at AAMI PB70 Level 3:

Impact Penetration	Hydrostatic Pressure	Weight	Fluid Resistance	Breathability
≤1.0 g	≥50 cm	Light		

Made with highly fluid resistant materials. Not PB70 rated.

SKU#	SIZE	GOWNS/BAG	BAGS/CASE	GOWNS/CASE
970B	UNIVERSAL	20	6	120
971B	EXTRA LARGE	20	6	120

① “Clean Release” Overhead Neck Closure

- Overhead neck closure promotes compliance and increases donning & doffing speed
- Removal of a shoulder seam creates a proper release at the back neck closure removing risk of contamination

② Soft & Strong Polyethylene Film

- Single, lightweight layer of protective polyethylene film provides a strong level of protection against fluids

③ Aligned Thumbloops

- Properly aligned thumbloops increases donning speed

