



LAB JACKET (HIP HIGH) KNIT CUFF • KNIT COLLAR

111XLJ / 112XLJ



Designed for comfort and mobility, VIVIT Supply Lab Jackets are constructed with breathable, lightweight SMS materials to keep wearers cool throughout the day. Its shorter, hip-length design allows for greater movement in clinical settings. Featuring a knit collar and cuffs, three versatile pockets to store accessories, and secure front snap closure, this jacket delivers practical protection for everyday wear.

SKU#	COLOR	SIZE	JACKETS/ BAG	BAGS/CASE	JACKETS/ CASE
1112LJ	BLUE	SMALL	10	3	30
1113LJ	BLUE	MEDIUM	10	3	30
1114LJ	BLUE	LARGE	10	3	30
1115LJ	BLUE	X-LARGE	10	3	30
1116LJ	BLUE	2X-LARGE	10	3	30
1122LJ	WHITE	SMALL	10	3	30
1123LJ	WHITE	MEDIUM	10	3	30
1124LJ	WHITE	LARGE	10	3	30
1125LJ	WHITE	X-LARGE	10	3	30
1126LJ	WHITE	2X-LARGE	10	3	30

① Functional Design

- Easy snap-front closures help secure the jacket while enabling quick donning and doffing
- Chest and waist pockets offer accessible storage for practical, everyday functionality

② Soft & Breathable SMS

- Specially designed material allows for the perfect mixture of fluid protection and comfort
- Maximizes breathability throughout the jacket

③ Extra Soft Knit Cuff

- Material provides maximum comfort during periods of extended use
- Snug fit to keep cuff in place over clothing



Not made with natural rubber latex



VIVIT Supply LLC
200 Continental Dr., Suite 401
Newark, DE 19713 USA

W vivitsupply.com

*Registered Trademark or Trademark of
PRIMED Medical Products Inc. or its affiliates.

V001