

Why Can't Tumor Markers Work Better?

You asked / We listened

VM-1 Tissue Marker

Future of Post Biopsy Care

A New Paradigm: From Detection to Lifelong Care

Screening has always been a cornerstone of breast cancer care—but long-term visibility is just as critical.

VM1 redefines what's possible after a biopsy, supporting ACR guidelines with multi-modality visibility and lifetime durability.



MRI Compatibility

The vial is filled with a specialized MRI contrast agent that significantly enhances visibility during scans. This agent maintains consistent brightness over time, facilitating accurate visualization of the marker within the body.

Non-metallic

Allergy free breast marker that delivers unmatched clarity and stability across all imaging modalities including MRI, Mammogram, Tomosynthesis, Ultrasound and PET.

Permanent Placement

Infused PEEK material with barium sulfate and carbon adheres securely to tissue, eliminating migration after implantation while producing a bright signal producing no artifacts.

The marker enhances visibility across all imaging modalities Ultrasound, Tomography, and Mammography for the life of the patient

