





HEATING WATER BATH

HB SERIES

Precision and performance for laboratories.

The new HB Series Water Baths are designed for high accuracy temperature adjustment. Ease-of use combined with precision simplifies workflows and maximize productivity.





ACCURATE AND DEPENDABLE

- With highly accurate temperature control systems HB Series Water Baths are designed to maintain water temperature on desired level from ambient to 100°C.
- Wide range of capacities from 5L to 25L*. The systems are ideal for all laboratory applications. With wide range of options, a suitable model is available for every usage.
- Simple maintenance for operators. Easily drain the water inside and clean the chamber with stainless steel interior.
- With safety features including audible alarm and over-heating protection, your water baths are long-lasting and dependable.
- Choose additional accessories* to your special needs and level-up your operations.

TECHNICAL DATA

	HB 305	HB 309	HB 313	HB 318
Power Input	1000 W	1000 W	1500 W	1500 W
Temperature Range	Room Temp-99°C	Room Temp-99°C	Room Temp-99°C	Room Temp-99°C
Temperature Control	Turning Knob	Turning Knob	Turning Knob	Turning Knob
Set Temperature Resolution	±1°C	±1°C	±1°C	±1°C
User Interface	Digital	Digital	Digital	Digital
Filling Volume Max.	5.5 L	9 L	13 L	18 L
Inner Height	65 mm	100 mm	150 mm	200 mm
Water Draining	Quick Coupling	Quick Coupling	Quick Coupling	Quick Coupling
Water Filling	On Top	On Top	On Top	On Top
Accessories	Standard*	Standard*	Standard*	Standard*
Material in Contact	Stainless Steel	Stainless Steel	Stainless Steel	Stainless Steel
Permissible Ambient Temperature	5-40°C	5-40°C	5-40°C	5-40°C
Permissible Relative Humidity	80%	80%	80%	80%
Voltage	220-240 V	220-240 V	220-240 V	220-240 V
Frequency	50/60 Hz	50/60 Hz	50/60 Hz	50/60 Hz

*Standard accessories include lid and standard tray. For additional accessories please contact for options.

**HB Series Water Baths can be designed according to your special needs (volume, dimensions, power input and accessories can be designed for your requirements). Contact to learn more about custom designed water baths.

FOR ALL YOUR QUESTIONS AND INQUIRIES

✉ support@debyetechnic.com

☎ (+90312) 485 40 42

☎ (+90312) 485 40 42

✉ info@debyetechnic.com

🌐 www.debyetechnic.com

in DEBYE TECHNIC



NOTE: Contents may be subject to modifications or corrections