



SHILEY EYE INSTITUTE
UC SAN DIEGO

Bio-interventional Cyclodialysis in Patients Undergoing Cataract Surgery and Failing IOP Control with Topical Medication Therapy


Robert N. Weinreb, MD, Leon W. Herndon, MD, Nathan M. Radcliffe,
MD, Sean A. Ianchulev, MD MPH

April 27, 2025


Paper 112558

Effective Interventions Needed for Moderate Glaucoma Patients with Medically Uncontrolled IOP

- Medical therapy for IOP reduction is highly effective, however, a subset of patients cannot be adequately controlled using medications alone
- Challenges with polypharmacy
Multi-medication adherence
Ocular surface
Continued progression
- 90% of MIGS used in mild-mod OAG in combination with phaco



Uveoscleral Outflow: A Validated Pathway





- **Uveoscleral outflow is an established therapeutic target**
- **Demonstrated Rx efficacy with pharmacotherapy (PGAs)**
- **Larger aqueous drainage capacity than trabecular pathway**
- **Negative oncotic gradient drives aqueous outflow**

Surgical Access to the Uveoscleral Outflow Pathway



Cyclodialysis History and Evolving Clinical Paradigm

1905	1920-1940+	1950	1960	2000+
First Cyclodialysis	Most common	Cyclodialysis	Supplanted by	Interventional ab-
Heine	glaucoma surgery	Ab-interno	trabeculectomy	interno approach ±
	External approach	approach		reinforcement



Cyclodialysis Evolution of Surgical Technique

Ab-externo >> Ab-Interno >> μ -interventional >> Reinforced



Material Matters!

Bio-tissue reinforcement of cyclodialysis

a novel bio-interventional approach for uveoscleral enhancement



Inventor
Dr. Ianchulev



Conforming Implant Material

Soft, scleral wall
compliant bio-
tissue; no vertical
rigid tip



No Hardware

Bio-conforming
soft tissue; no
plastic, metal or
rigid foreign body




Goniometric Controlled Depth Implantation

Transparent
goniometric tip to
preempt
anteriorized
deployment



No Rebound Movement


Enhanced post-
deployment
fixation - tissue
expansion in cleft



Blakney AK, Swartzlander MD, Bryant SJ. The effects of substrate stiffness on the in vitro activation of macrophages and in vivo host response to poly(ethylene glycol)-based hydrogels. J Biomed Mater Res A. 2012 Jun;100(6):1375-86. doi: 10.1002/jbm.a.34104. Epub 2012 Mar 7. PMID: 22407522; PMCID: PMC3339197.

Material Matters!

100% Biologic, hardware-free, allogeneic reinforcement of cyclodialysis





Homologous Allograft

Minimal foreign body reaction
Designed for corneal safety

Porous

For aqueous conductivity



Foreign body response: Reduced Macrophage activation with more homologous implant material biomechanics

	Tensile Strength MPa	SC Implant vs Tissue Mismatch
Scleral Allograft	1-2 ¹	0
Titanium	240 ³	239
Polyimide	80 ²	70

	Cornea Permeability (cm/s)	Sclera Permeability (cm/s)
BENZOLAMIDE	1.4 E-7	2.0E-5
INULIN	5.5 E-7	9.0 E-6
PROPRANOLOL	3.1 E-6	5.8 E-5
SUCROSE	4.4 E-6	2.2 E-5

Micro-interventional Cyclodialysis with Scleral Reinforcement




Study Design and Baseline Characteristics

- CREST is a multicenter, interventional, safety and effectiveness, real world clinical study (n=410)
- Medically uncontrolled cohort of glaucoma patients with a confirmed diagnosis of open-angle glaucoma, operable cataract, baseline IOP ≥ 21 mmHg. Prospectively followed for 24 months
- Interventions: interventional ab-interno cyclodialysis with subsequent bio-tissue scleral reinforcement with or without phaco surgery
- Analysis: Ocular safety, concomitant procedures, IOP and medication usage over 12 months

Baseline Characteristics	
Sample size, eyes, N	51
Patients	38
Age, mean \pm SD, years	71.4 \pm 10.5
Gender: female, %	44.7%
Pseudophakic	68%
Baseline BCVA, mean \pm SD	0.40 \pm 0.32
Baseline IOP, mmHg, mean \pm SD	25.7 \pm 4.4
Number of IOP-lowering drugs, mean \pm SD	1.2 \pm 1.3

Efficacy Outcomes



	BL (n=51)	1M (n=50)	6M (n=41)	12M (n=45)
IOP, mean	25.7	15.7	16.0	15.2
IOP, SD	4.4	4.4	4.0	4.5

	BL (n=51)	1M (n=50)	6M (n=41)	12M (n=45)
Mean IOP Meds	1.2	0.9	0.9	0.7
+SD	1.3	1.2	0.9	0.9

Safety Events

	N=51
BCVA, decimal, mean \pm SD	0.73 \pm 0.23
IOP increase (10 mmHg or >30mmHg), n (%)	4 (7.8)
IOP increase (>30mmHg), n (%)	4 (7.8)
>2 lines drop in BCVA, n (%)	0 (0)
Significant intraoperative complications (bleeding, tissue damage)	0 (0)
Persistent inflammation (>1M), n (%)	0 (0)
Severe inflammation (grade 4+), n (%)	0 (0)
Persistent hyphema (>1M), n (%)	1 ¹ (1.9)
Severe post-operative hyphema (>3 mm), n (%)	0 (0)
Persistent corneal edema (>1M), n (%)	0 (0)
Bio-tissue migration, n (%)	0 (0)
Cystoid macular edema	2 ² (3.9)
IOP <6 mm on more than 2 visits, without maculopathy	0 (0)
Hypotony with Maculopathy	0 (0)
Post-operative Laser Enhancement	2 ³ (3.9)
Secondary Glaucoma Intervention	4 ³ (7.8)

1. Fully resolved
2. With phaco surgery and diabetic retinopathy. No hypotony
3. 1 case of aqueous humor release (wound bupring) at the 1-week postop visit, 1 case of Ahmed tube, 1 Case of Xen implant, and 1 case of CPC/Goniotomy

Conclusion

- Bio-interventional enhancement of uveoscleral outflow is achieved through ab-interno cyclodialysis and homologous allograft scleral reinforcement
- Allograft tissue demonstrates inert, biocompatible, and bio-conforming material properties
- Allograft bioscaffold shows durable structural stability with no post-operative migration.
- Bio-interventional cyclodialysis surgery demonstrates robust IOP lowering effect in the setting of cataract surgery in medically uncontrolled OAG patients

