



myoncare

Exporting Answers

Version 1.12





PDF Report



Export symbol: patient responses to Caretask components can be exported in CSV or PDF format for a selectable date range.

myoncare Careplan Manager allows healthcare professionals to decide if their patients get a PDF report of conducted Caretasks. Through the report, the patients can review the results and forward the information to another doctor to assess it, if needed.

The PDF Report contains a generic first page which includes general information about the patient like name, Patient-ID, date of birth, ICD-Code, gender, the institution name and address of the institution. Also, it includes the name of the Caretask, the description of the Caretask (if available), the date when the report was generated, total recurrences and lifespan of the Caretask.

The following pages of the report include a header with the name of the Caretask and when the report was generated. Besides that, it always contains the components and the given answers. The components include an index number and the title, as well as all answer options (depending on the type of component).


The PDF reports include both section scores and overall scores, providing deeper insights into specific areas of patient progress. This ensures that healthcare professionals and patients have access to more comprehensive and actionable data. By offering detailed section-level analysis, the reports support better care coordination, more informed decision-making, and improved communication between patients and their healthcare providers.


Note: Only the components that were displayed to the patient while answering the caretask are available in the report.

Components which include a defined jump for one or more options will include an indicator which displays which option leads to which of the following components.

To generate the PDF-Report, the healthcare professionals have to activate the option 'Create a Caretask report for the patient' within the Caretask settings.

Once activated, additional settings will open. The healthcare professional has to decide, if the report should be generated every time a patient or healthcare professional answers a caretask or only if a threshold is exceeded. Besides that, the healthcare professional has to decide in which language the Caretask report should be generated. Per default the language is set in the same way as the healthcare professional has set the myoncare Careplan Manager. The language setting will only affect the generic information of the report, while the content will not be translated.





Additionally, healthcare professionals can review the report by enabling the "Report review" toggle. They also have the option to select a 'Default Review Caretask' that will automatically be associated with the reports generated. These reports can be easily accessed by clicking on the "Ready-to-Review Reports" notification icon located at the top right of the screen.

The patient gets a push notification, through the PatientApp, when a report is generated. The patient can access the report either clicking on the push notification, or going to the file section, within the Mobile. The generated report will be also found within the file tab of the patient profile in the myoncare Careplan Manager.

Note: Reports that require review are not available in the file sections and are not visible to the patients.

The patient gets a push notification, through the Mobile, when a report is generated. The patient can access the report either clicking on the push notification, or going to the file section, within the Mobile. The generated report will be also found within the 'File' tab of the patient profile in the myoncare Careplan Manager.

Note: Reports that require review are not available in the 'File' section and are not visible to the patients.



We make healthcare scalable.

www.myon.academy | www.myoncare.com

Bei Fragen: service@myoncare.com