



myoncare

Creating a Patient

Version 1.12





Creating a patient

NOTE: In general, if a myoncare Careplan Manager installation includes an HL7 interface (non-mandatory module), parts of the patient data stored in the hospital information system (HIS) can automatically be imported into myoncare Careplan Manager and do not have to be entered manually.

For creating a patient manually within the myoncare Careplan Manager, the healthcare professional needs to do the following steps:



- Click on the action button (“+”-button) within the patient overview.



- Click on the “Create patient”-button.
- Add the patient information in the mask accordingly. All input fields which are marked with a “*” are mandatory (First name, Last name, and Date of birth).

NOTE: Patient-ID is the unique identifying attribute to other systems (HIS/KIS) and needs to be identical to the medical record within interfacing systems. Reference-ID is an open pseudonym that each site can choose for its patients. It is not necessarily unique.

Medical information such as patient ID, case number, OPS Reference or reference ID can be adjusted in the myoncare Careplan Manager at any time.

In case no specific Patient-ID is defined for a patient within the creation process, the patient receives a standardized combination of 3 characters of the first name, last name and date of birth.



We make healthcare scalable.

www.myon.academy | www.myoncare.com

Bei Fragen: service@myoncare.com