


TiLINK™

SI FUSION SYSTEM


 SURGENTEC®



TiLINK™

SI FUSION SYSTEM

TiLINK-L SI Fusion System is intuitively engineered to help compress the sacroiliac joint and provide an optimal platform for fusion. TiLINK-L has capabilities to accomplish both lateral or posterior oblique approaches.



Nanotex® Surface Technology

The hydrophilic surface created by Nanotex® takes advantage of ionic bonds to attract and hold BMA & other fluids that may contain key growth factors

Bone Ingrowth Windows

Implant contains strategically placed windows for potential points of bony ingrowth

Harvesting Channel

The channel gathers autograft as the screw is inserted to help facilitate fusion

Self-Tapping Tip

Self-tapping tip forms threads in bone to create a path and secure the implant

Cannulated Screw

Cannulated screw allows for precise placement of implant

SCREW SIZES:

All Screws are 10mm in diameter.

