

Order at: orders@nonin.com or fax: +1 763-577-5521

Customer Service inquiries: customerservice@nonin.com or (800) 356-8874



MODEL NUMBER	DESCRIPTION	CABLE LENGTH	APPLICATION	PACKAGING	WristOx ₂ [®] Model 3150	CO-Pilot [™] **	PalmSAT [®] 2500/2500A	8500 Series	Model 7500**	Avant [®] 9700
Reusable Soft Sensors 8000S Series										
 8000SS 8000SS-3	Small Soft Sensor 7.5 – 12.5 mm digit 0.3 – 0.5 inches digit	1m 3m	Index, middle, or ring finger	Single			✓	✓	✓	✓
8000SS-WO2	Small Soft Sensor – WristOx₂ 7.5 – 12.5 mm digit 0.3 – 0.5 inches digit	0.3m	Index, middle, or ring finger	Single	✓					
 8000SM 8000SM-3	Medium Soft Sensor 10 – 19 mm digit 0.4 – 0.75 inches digit	1m 3m	Index, middle, or ring finger	Single			✓	✓	✓	✓
8000SM-WO2	Medium Soft Sensor – WristOx₂ 10 – 19 mm digit 0.4 – 0.75 inches digit	0.3m	Index, middle, or ring finger	Single	✓					
 8000SL 8000SL-3	Large Soft Sensor 12.5 – 25.5 mm digit 0.5 – 1.0 inches digit	1m 3m	Index, middle, or ring finger	Single			✓	✓	✓	✓
8000SL-WO2	Large Soft Sensor – WristOx₂ 12.5 – 25.5 mm digit 0.5 – 1.0 inches digit	0.3m	Index, middle, or ring finger	Single	✓					
Reusable Finger Clip, Ear Clip and Reflective Sensors 8000 Series										
 8330AA	Finger Clip Sensor – Adult > 30 kg / > 66 lbs	1m	Index, middle, or ring finger	Single		✓				
 8000AA 8000AA-2 8000AA-3	Finger Clip Sensor – Adult > 30 kg / > 66 lbs	1m 2m 3m	Index, middle, or ring finger	Single			✓	✓	✓	✓
8000AA-WO2	Finger Clip Sensor – WristOx₂ > 30 kg / > 66 lbs	0.3m	Index, middle, or ring finger	Single	✓					
 8000AP 8000AP-3	Finger Clip Sensor – Pediatric 8 – 30 kg / 18 – 66 lbs	1m 3m	Index, middle, or ring finger	Single			✓	✓	✓	✓
 8000Q2	Ear Clip Sensor > 40 kg / > 88 lbs	1m	Earlobe	Single	✓		✓	✓	✓	✓
8000R 8000H	Reflectance Sensor > 30 kg / > 66 lbs Sensor Holder 10 caps / 20 stickers	1m N/A	Forehead	Single			✓	✓	✓	✓
Reusable Flex Sensors with FlexiWrap 8000 Series										
 8000J 8000J-3 8000JFW	Reusable FlexSensor – Adult > 22 kg / > 44 lbs incl. 25 Adult FlexiWraps	1m 3m N/A	Index, middle, or ring finger	Single 24 per pack			✓	✓	✓	✓
8000J-WO2 8000JFW	Reusable FlexSensor – WristOx₂ > 22 kg / > 44 lbs incl. 25 Adult FlexiWraps	0.3m N/A	Index, middle, or ring finger	Single 24 per pack	✓					
 8008J 8008JFW	Reusable FlexSensor – Infant 2 – 22 kg / 4.4 – 44 lbs incl. 25 Infant FlexiWraps	1m N/A	Toe	Single 24 per pack			✓	✓	✓	✓
 8001J 8001JFW	Reusable FlexSensor – Neonatal > 2 kg / 4.4 lbs incl. 25 Neonatal FlexiWraps	1m N/A	Side of foot	Single 24 per pack			✓	✓	✓	✓
Single-Patient-Use Cloth Sensors 6000C Series										
 6000CA	Cloth Sensor – Adult > 30 kg / > 66 lbs	1m	Index, middle, or ring finger	Single 24 per pack			✓	✓	✓	✓
6000CA-WO2	Cloth Sensor – Adult > 30 kg / > 66 lbs	0.3m	Index, middle, or ring finger	Single 15 per pack	✓					
 6000CP	Cloth Sensor – Pediatric 10 – 40 kg / 22 – 88 lbs	1m	Index, middle, or ring finger	Single 24 per pack			✓	✓	✓	✓
 6000CI	Cloth Sensor – Infant 2 – 20 kg / 4.4 – 44 lbs	1m	Big toe	Single 24 per pack			✓	✓	✓	✓
 6000CN	Cloth Sensor – Neonatal/Adult Neonatal: < 2 kg / < 4.4 lbs Adult: > 30 kg / > 66 lbs	1m	Foot	Single 24 per pack			✓	✓	✓	✓
Single-Patient-Use Flexi-Form[®] III Sensors 7000 Series										
 7000A	Flexi-Form Sensor – Adult > 30 kg / > 66 lbs	1m	Index, middle, or ring finger	Single 24 per pack			✓	✓	✓	✓
 7000P	Flexi-Form Sensor – Pediatric 10 – 40 kg / 22 – 88 lbs	1m	Index, middle, or ring finger	Single 24 per pack			✓	✓	✓	✓
 7000I	Flexi-Form Sensor – Infant 2 – 20 kg / 4.4 – 44 lbs	1m	Big toe	Single 24 per pack			✓	✓	✓	✓
 7000N	Flexi-Form Sensor – Neonatal/Adult Neonatal: < 2 kg / < 4.4 lbs Adult: > 30 kg / > 66 lbs	1m	Foot	Single 24 per pack			✓	✓	✓	✓
Accessories										
 UNI-EXT-1 UNI-EXT-3 UNI-EXT-6 UNI-EXT-9	Extension Cable	1m 3m 6m 9m		Single			✓	✓	✓	✓
 UNI-RA-3 UNI-RA-0	Right Angle Connector	3m 15cm		Single			✓	✓	✓	✓
 3150I	Adapter Cable to 9-pin (D-Connector). Connects PureLight [®] sensors to WristOx ₂ Model 3150	12cm		Single	✓					

Authorized Distributor:

Nonin Medical, Inc.

13700 1st Avenue North
Plymouth, MN • 55441-5443 • U.S.A.

Tel: +1.763.553.9968 1.800.356.8874

Fax: +1.763.577.5521

Email: info@nonin.com

Nonin Medical B.V.

Doctor Paul Janssenweg 150
5026 RH Tilburg • Netherlands

Tel: +31 (0)13-45 87 130

Email: infointl@nonin.com