



**ApneaLink™ Air**  
Home sleep testing device



**ApneaLink™ Air**  
Home sleep testing device

## ApneaLink Air complete set

Includes oximetry accessories and EasySense respiratory effort sensor.



Oximetry Set with Reusable Soft Sensor	1 XPOD LP Oximeter 1 XPOD LP Oximeter Clip 1 NONIN Oximeter Reusable Soft Sensor (8000SM)	<b>22370*</b>
Oximetry Set with Disposable Soft Sensor	1 XPOD LP Oximeter 1 XPOD LP Oximeter Clip 3 NONIN Oximeter Single-use Soft Sensors (7000A)	<b>22380*</b>
Nasal Cannulas	Nasal Cannulas – 25 pk	<b>70388</b>
Other accessories	ApneaLink Air Reusable Belt ApneaLink Air Carry Bag ApneaLink Air Software CD 10.1 USB Cable	<b>629052</b> <b>22373</b> <b>22372</b> <b>22375</b>

### Technical specifications

Recorder	Enhanced Hardware EasySense respiratory effort sensor Supports oximetry Light indicators - Test complete indicator - Respiratory flow indicator - Effort Sensor connection indicator - Oximetry connection indicator
Signal recording	Respiratory effort Respiratory flow Snore Blood oxygen saturation Pulse Battery voltage
Download application	PC download
Internal battery	2 AAA alkaline batteries
Internal memory	Recording period: 48 hours Internal clock
Dimensions	Recorder: 60 x 86 x 28 mm (2.4" x 3.4" x 1.4") Pulse oximeter: 53 x 20 x 15 mm (2.1" x 0.8" x 0.6") Recorder weight: 66 g (2.3 oz)

Please contact your ResMed representative for more details.

### Ordering information and product codes

ApneaLink Air Complete Set	ApneaLink Air recorder device Respiratory effort sensor Software installation CD USB download cable 1 reusable belt 3 nasal cannulas	ApneaLink Air patient instructions XPOD LP Oximeter XPOD LP Oximeter clip 1 Nonin oximeter reusable soft sensor 8000SM Carry bag 2 AAA batteries	<b>22354</b>
----------------------------	---	---	--------------

### Accessories & consumables

Oximetry components	XPOD LP Oximeter	<b>22374*</b>
	XPOD LP Oximeter Clip	<b>22371*</b>
	NONIN 7000A single use oximetry sensors - 24 pk	<b>22338</b>
	NONIN 6000A single use oximetry sensors - 24 pk	<b>22337</b>
	NONIN Oximeter Reusable Soft Sensor (8000SS) - Small	<b>70567</b>
	NONIN Oximeter Reusable Soft Sensor (8000SM) - Medium	<b>70568</b>
	NONIN Oximeter Reusable Soft Sensor (8000SL) - Large	<b>70413</b>
	NONIN Oximeter Reusable Finger Clip Sensor - 8000AA (1m)	<b>631803</b>

\* Available from June 2014 onwards



## Home sleep testing made easy

ResMed's new ApneaLink™ Air provides healthcare professionals with a portable and easy-to-use home sleep testing device for patients who might have sleep disordered breathing. ApneaLink Air is a cost effective and convenient solution to test patients for sleep apnea in the comfort of their own bedroom. Testing patients in their own environment can help improve their experience and acceptance of the world of sleep apnea therapy right from the start.



**DISTRIBUTED BY**  
ResMed Pty Ltd  
1 Elizabeth Macarthur Drive Bella Vista NSW 2153 Australia

ApneaLink is a trademark of ResMed R&D Germany GmbH and is registered in the U.S. Patent and Trademarks Office. See [www.resmed.com](http://www.resmed.com) for other ResMed locations worldwide. Specifications may change without notice. © 2014 ResMed. 1017396/1 2014-04

[ResMed.com](http://ResMed.com)

[ResMed.com](http://ResMed.com)

