

Mediclose

Close 12mm trocar ports easily and efficiently with the Mediclose Trocar Port Closure System. Developed by a Surgeon for Surgeons, Mediclose reduces the time and effort to close, by moving the fascia into the path of the superficial-facing needles. With suture and needle self-contained in each cartridge, the spreading mechanism elongates the wound, so the surface-directed needle deployment gets a consistent and organ-sparing bite.

SPECIFICATION

Trocar Port Size: 12mm

Delivery Mechanism: Spring-loaded

Suture & Needles: All-inclusive

KEY FEATURES

- Suturing time <30 seconds
- One-handed use
- Superficial needle orientation
- Quick-change cartridge



COMPLIANCE & STANDARDS

- FDA: Approved

Our Office:

300 W Antelope Road
White City, OR 97503

Contact Information:

Tel: 541-770-5905
Mail: info@lwshealthandsafety.com