

# Foundation Model Distillation to Develop Robust and Efficient AI Models for MSI Testing in Endometrial Cancer



September 2025

Alexandre Filiot<sup>1</sup>, Elin Samuelsson<sup>1</sup>, Mélanie Lubrano<sup>1</sup>, John Klein<sup>1</sup>, Aurore Lyon<sup>1</sup>, Corentin Wiscart<sup>1</sup>, Arnaud Fouillet<sup>1</sup>, Lydwine van Praet<sup>1</sup>, Aymeric Vinçotte<sup>1</sup>, Catherine Genestie<sup>2</sup>, Frederik Deman<sup>3</sup>, Luis Cano Ayestas<sup>1</sup>, Katharina von Loga<sup>1</sup>

1 Owkin, Paris, France, 2 Gustave Roussy, Department of Pathology, Paris, France, 3 Ziekenhuis aan de stroom (ZAS), Department of Pathology, Antwerp, Belgium

## Background

→ AI models for predicting microsatellite instability (MSI) from H&E slides require diverse, multi-center datasets. This is challenging in cancers with lower MSI prevalence, such as endometrial cancer (EC). While recent foundation models (FM) address data scarcity, they often fail to overcome biases in digital pathology and contribute to high computational costs, limiting clinical applicability.


## Objectives

→ This study demonstrates that a distilled version of a large foundation model can offer an efficient, robust AI-assisted diagnostic tool for MSI testing in EC.

## Method

→ **Task:** Predict MSI status from digitized H&E slides in GC and EC

→ **Model:** Feature extraction with Foundation model + MIL (Chowder)



→ **Foundation model comparison:** We compared the performance of a large foundation model ("H0") [1] with a 13x-smaller version ("H0-mini"). [2]

→ **Data:** multi-center cohorts

Indication	Cohort	N. patients (MSI-h)	N WSI	Tissue type
EC	TCGA-UCEC	478 (153)	559	Resections
EC	ZNA-GZA	130 (41)	216	Resections + biopsies
EC	GR	109 (41)	109	Resections

## Discussion

→ Advanced distillation techniques yield an enhanced, efficient AI model for MSI detection in EC, improving detection capabilities while lowering computational costs. These findings highlight the potential for broader implementation of robust MSI testing in routine clinical practice, aiding streamline diagnosis and targeted therapy in this under-researched cancer population.

## References & acknowledgements

[1] Saillard, C., Jenatton, R., Llinares-López, F., Mariet, Z., Cahané, D., Durand, E., Vert, J.P.: H-optimus-0 (2024), <https://github.com/biopotimus/releases/tree/main/models/h-optimus/v0>

[2] Filiot, A., Dop, N., Tchita, O., Riou, A., Dubois, R., Peeters, T., ... & Olivier, A. (2025). Distilling foundation models for robust and efficient models in digital pathology. arXiv preprint arXiv:2501.16239.

The results published here are in whole or part based upon data generated by the TCGA Research Network. This work was granted access to the HPC resources of IDRIS under the allocation AD011012519 made by GENCI. This work was funded as part of a collaboration between Owkin and MSD.

## Results

→ **Performance:** The H0-mini model achieved an AUROC of 0.769 (0.738, 0.760, 0.810) and specificity of 0.497 (0.451, 0.393, 0.648) at 0.9 sensitivity across three left-out cohorts, outperforming H0 by 0.027 ( $p=0.03$ ) and 0.291 ( $p<0.001$ ), respectively.

Evaluation	Training on TCGA-UCEC		Training on GR and ZNA-GZA
	GR	ZNA-GZA	TCGA-UCEC
H0	AUC: 0.810 (0.733, 0.883) Sp: 0.493 (0.372, 0.609)	AUC: 0.760 (0.691, 0.827) Sp: 0.300 (0.231, 0.375)	AUC: 0.738 (0.696, 0.780) Sp: 0.363 (0.314, 0.412)
Distilled	AUC: <b>0.842</b> (0.767, 0.912) Sp: <b>0.648</b> (0.526, 0.754)	AUC: <b>0.768</b> (0.704, 0.832) Sp: <b>0.393</b> (0.316, 0.475)	AUC: <b>0.761</b> (0.721, 0.803) Sp: <b>0.451</b> (0.404, 0.501)

**Table: MSI classification metrics on external cohorts for H0 and H0-mini**

→ **Robustness:** H0-mini achieved an intra-class correlation of 0.725 (95% CI: 0.651-0.793) between paired resections and biopsies, compared to 0.665 (95% CI: 0.573-0.750) for H0 ( $p<0.01$ ).

Bootstrap over 1000 repeats	Pearson R	CCC (Concordance Correlation Coefficient)	ICC (Intra-Class Correlation)
H0	0.668 (0.540, 0.775)	0.614 (0.482, 0.726)	0.665 (0.573, 0.750)
Distilled	<b>0.727 (0.625, 0.810)</b>	<b>0.718 (0.612, 0.802)</b>	<b>0.725 (0.651, 0.793)</b>
<b>t-test p-value from paired bootstrap distributions</b>			
	<b>0.043</b>	<b>0.003</b>	<b>0.009</b>

**Table: Concordance between Rx and Bx scores for endo lock model depending on the FM**

→ **Computational cost:** The H0-mini model is 9.1 times faster than H0 on a GPU machine. On standard pathology workstation, processing a WSI takes 15.2 minutes with H0-mini, compared to 120.4 minutes for H0, highlighting the intractability of large FMs without dedicated infrastructure.

Tiling time in minutes	H0 GPU	Distilled GPU	H0 CPU	Distilled CPU	Ratio GPU averaged over slide	Ratio CPU	Ratio Distilled CPU / H0 GPU
10 TCGA slides (152k tiles in total)	84	9	-	-	9.0 (0.8)	-	-
2 TCGA slides (29.5k tiles in total)	-	-	227	35	-	6.7 (0.4)	1.9 (0.5)

**Table: Tiling time in minutes on GPU and CPU**