


# Biodesign®

OTOLOGIC REPAIR GRAFT

Harvest *results*,  
not patient tissue<sup>1</sup>



## Reliable Closure

Biodesign® material remodels into natural host tissue with an overall success rate of 91% across published literature<sup>1-9</sup> with no statistically significant difference in audiometric results when compared to temporalis fascia.<sup>1,10</sup>

## Excellent Handling

Biodesign® material is easy to manipulate, allowing for improved surgical precision during graft placement.<sup>1</sup>


## Time Saving

The Biodesign® Otolgic Repair Graft reduces the need to harvest autologous tissue, significantly decreasing intraoperative time.<sup>1</sup>

## Biodesign® Otologic Repair Graft

### Available product sizes

Shown at actual size.



Order Number	Reference Part Number	Size cm	Nominal Thickness mm
G44840	ENT-OTO-0.4-0.6-2	0.4, 0.6	0.25
G44839	ENT-OTO-0.6-0.9-2	0.6, 0.9	0.25
G44451	ENT-OTO-2.5X2.5-2	2.5 x 2.5	0.25
G44452	ENT-OTO-5X5-2	5.0 x 5.0	0.25

### References

- D'Eredità R. Porcine small intestinal submucosa (SIS) myringoplasty in children: a randomized controlled study. *Int J Pediatr Otorhinolaryngol.* 2015;79(7):1085-1089.
- Cass ND, Hebbe AL, Meier MR, et al. Pediatric primary tympanoplasty outcomes with autologous and non-autologous grafts. *Otol Neurotol.* 2022;43(1):94-100.
- Chen CK, Hsieh LC. Clinical outcome of exclusive endoscopic tympanoplasty with porcine small intestine submucosa in 72 patients. *Clin Otolaryngol.* 2020;45(6):938-943.
- Barron C, Lukens J, Niermeyer W, et al. Investigation of novel grafts in use for pediatric tympanoplasty. *Ann Otol Rhinol Laryngol.* 2019;128(12):1111-1115.
- Redaelli De Zinis LO, Berlucchi M, Nassif N. Double-handed endoscopic myringoplasty with a holding system in children: preliminary observations. *Int J Pediatr Otorhinolaryngol.* 2017;96:127-130.
- James AL. Endoscope or microscope-guided pediatric tympanoplasty? Comparison of grafting technique and outcome. *Laryngoscope.* 2017;127(11):2659-2664.
- Ranguis SC, Leonard CG, James AL. Prospective comparison of pediatric endoscopic lateral graft and interlay tympanoplasty. *Otol Neurotol.* 2021;42(6):867-875.
- Wang N, Isaacson G. Collagen matrix as a replacement for Gelfilm for post-tympanostomy tube myringoplasty. *Int J Pediatr Otorhinolaryngol.* 2020;135:110136.
- Yawn RJ, Dedmon MM, O'Connell BP, et al. Tympanic membrane perforation repair using porcine small intestinal submucosal grafting. *Otol Neurotol.* 2018;39(5):e332-e335.
- Dontu P, Shaigany K, Eisenman DJ. Anatomic and audiometric outcomes of porcine intestinal submucosa compared to autologous fascia for tympanic membrane repair. Poster presented at: Combined Otolaryngology Spring Meetings, COSM 2022; April 27 – May 1, 2022; Dallas, TX.

® indicates U.S. trademark registration. All trademarks and/or images are the property of their respective owners or holders.



## Biodesign®

OTOLOGIC REPAIR GRAFT

For more detailed product information, including the Instructions for Use, which contains the indication statement, contraindications, precautions, and potential complications, visit [www.everismedical.com](http://www.everismedical.com).

