

COMPATIBILITY GUIDE




-  **Fully assembled** Mass Spectrometer Bench.
-  **Fully movable** with lockable castor wheels.
-  Integrated vacuum pump(s) **noise enclosure**.
-  Fan **heat extraction** system.
-  Hoses & wires **grommets**.
-  Vacuum pump(s) **dampening vibration** system.
-  Overheating **temperature alarm**.




Delivered fully assembled




Fully movable



Integrated noise enclosure



Fan heat extraction system
Hoses & wires grommets



Dampening vibration system



Overheating temperature alarm

COMPATIBILITY GUIDE

TECHNOLOGY	WATERS INSTRUMENTS	ionBench MODELS
 MS only	BioAccord, SQD-2, TQD, TQ-S Micro, G2-XS Q-ToF	BCH120NE58
	TQ-S	BCH120NE78
	Xevo G3 Q-ToF	BCH13097-G3
Recommended by Waters site prep guide	Xevo MRT & CPU Holder for EPC included	BCH87123-MRT



TQD
BCH120NE58



Xevo MRT & LC Acquity
BCH87123-MRT & BCHLC4575 suggested



Xevo G3 Q-ToF
BCH13097-G3

COMPATIBILITY GUIDE

TECHNOLOGY	WATERS INSTRUMENTS	ionBench MODELS
LC/MS 	TQ-S Micro, TQD, SQD-D & LC Acquity - Short version	BCH120NE58
	TQ-XS & LC Acquity - Short version	BCH13097-TQXS
	TQ-Absolute & LC Acquity - Short version	BCH13097-TQA
	BioAccord, SQD-2, TQD, TQ-S Micro, G2-XS Q-ToF	BCH190NE58
	TQ-S & LC Acquity	BCH190NE78
	TQ-XS & LC Acquity - Wide version	BCH19097-TQXS
	TQ-Absolute & LC Acquity - Wide version	BCH19097-TQA
	Xevo G3 Q-ToF & LC Acquity	BCH19097-G3



BioAccord & LC Acquity
Short version
BCH120NE58



TQ-S & LC Acquity
BCH190NE78



TQ-Absolute & LC Acquity
BCH19097-TQA

COMPATIBILITY GUIDE



 **Fully assembled** Liquid Chromatography Bench.

 **Fully movable** with lockable castor wheels.


↑↓ Electrically **height adjustable**.

 **Chemical resistant** worksurface.

 **Tubing management**.

 **Minimize dead volume** between LC and MS.

 **Safe access** to mobile phase bottles.




Delivered fully assembled



Fully movable



Weight capacity up to 250 kg



Tubing path system
for waste management



Digital display for height
adjustment



Chemically resistant and easy to
clean.

COMPATIBILITY GUIDE

TECHNOLOGY	WATERS INSTRUMENTS	ionBench MODELS
LC only	LC Acquity - 1 Stack	BCHLC4575
	Alliance iS	BCHLC6765
	LC Acquity - 2 Stacks, Acquity & Sample Organizer	BCHLC9075



LC Acquity
BCHLC4575



**LC Acquity &
Sample Organizer**
BCHLC9075



Alliance iS
BCHLC6765

CONTACT US



EUROPE

+33 (0) 9 70 44 65 57

contact@ionbench.com



US

+1 984-317-9236

contact_us@ionbench.com