

COMPATIBILITY GUIDE




-  **Fully assembled** Mass Spectrometer Bench.
-  **Fully movable** with lockable castor wheels.
-  Integrated vacuum pump(s) **noise enclosure**.
-  Fan **heat extraction** system.
-  Hoses & wires **grommets**.
-  Vacuum pump(s) **dampening vibration** system.
-  Overheating **temperature alarm**.




Delivered fully assembled




Fully movable



Integrated noise enclosure



Fan heat extraction system
Hoses & wires grommets



Dampening vibration system



Overheating temperature alarm

COMPATIBILITY GUIDE

TECHNOLOGY	SCIEX INSTRUMENTS	ionBench MODELS
MS only	Triple Quad 3500, QTRAP 4500, Triple Quad 4500, Triple Quad 5500+ QTRAP Ready, X500B, X500R	BCH120NE58
	QTRAP 6500+, Triple Quad 6500+	BCH120NE78
	Triple Quad 7500 QTRAP Ready	BCH110100-SCIEX7500
	ZenoTOF 7600	BCH16097-ZENOTOF-7600
	Echo MS+	BCHLC11070_SCIEXSET p/n 5321658
LC/MS	Triple Quad 3500, QTRAP 4500, Triple Quad 4500, Triple Quad 5500+ QTRAP Ready, X500B, X500R & ExionLC Series	BCH190NE58
	QTRAP 6500+, Triple Quad 6500+ & ExionLC Series	BCH190NE78
	Triple Quad 7500 QTRAP Ready & ExionLC Series	BCH19097-SCIEX7500



**Triple Quad 7500 QTRAP Ready
&
ExionLC Series**
BCH19097-SCIEX7500



ExionLC & ZenoTOF 7600
BCHLC6765 & BCH190NE58-LEFT



Echo MS+
BCHLC11070_SCIEXSET
p/n 5321658

COMPATIBILITY GUIDE

TECHNOLOGY	SCIEX INSTRUMENTS	ionBench MODELS
MS & Dry pump	Triple Quad 7500 QTRAP Ready	BCH110100-SCIEX7500-DRY
	QTRAP 6500+, Triple Quad 6500+	BCH13097-SCIEX6500-DRY
	ZenoTOF 8600 system	BCH17097-ZENOTOF-8600
	ZenoTOF 7600	BCH16097-ZENOTOF-7600-DRY
LC/MS & Dry pump	QTRAP 6500+, Triple Quad 6500+ & ExionLC Series	BCH19097-SCIEX6500-DRY
	Triple Quad 7500 QTRAP Ready & ExionLC Series	BCH19097-SCIEX7500-DRY



Triple Quad 6500 & Dry pump
BCH13097-SCIEX6500-DRY



**ExionLC & Triple Quad 7500 QTRAP
Ready & Dry pumps**
BCH19097-SCIEX7500-DRY




Triple Quad 7500 QTRAP
BCH13097-SCIEX7500-DRY

COMPATIBILITY GUIDE



-  **Fully assembled** Liquid Chromatography Bench.
-  **Fully movable** with lockable castor wheels.
-  Electrically **height adjustable**.
-  **Chemical resistant** worksurface.
-  **Tubing management**.
-  **Minimize dead volume** between LC and MS.
-  **Safe access** to mobile phase bottles.




Delivered fully assembled



Fully movable



Weight capacity up to 250 kg



Tubing path system
for waste management




Digital display for height
adjustment



Chemically resistant and easy to
clean.

COMPATIBILITY GUIDE

TECHNOLOGY	SCIEX INSTRUMENTS	ionBench MODELS
LC only	ExionLC, ExionLC 2.0, ExionLC 2.0+	BCHLC6765
	ExionLC AD	BCHLC9075



ExionLC
BCHLC6765



ExionLC AD
BCHLC9075

CONTACT US



EUROPE

+33 (0) 9 70 44 65 57

contact@ionbench.com



US

+1 984-317-9236

contact_us@ionbench.com