

Travelling with Wheelchairs and Mobility Aids

To ensure your safety and comfort, special rules apply when travelling with battery-powered wheelchairs and mobility aids. Please review the information below before your journey.

Wheelchairs with Spillable Batteries

- Wheelchairs or mobility aids with spillable (liquid) batteries are not accepted for carriage.

Wheelchairs with Lithium Batteries

Wheelchairs and mobility aids powered by lithium batteries may be accepted under the following conditions:

- The battery may:
 - Remain installed, if it is securely attached and protected by the design of the device, or
 - Be removed, following the manufacturer's instructions
- If the battery is removed:
 - It must be carried in the cabin (hand baggage only)
 - Maximum allowed:
 - 1 battery up to 300 Wh, or
 - 2 batteries up to 160 Wh each
- Spare batteries:
 - Allowed up to a total of 300 Wh per mobility aid
 - Must be individually protected (e.g. in protective pouches)
 - Must be carried in the passenger cabin only

Wheelchairs with Non-spillable wet, nickel-metal hydride, or dry batteries

These mobility aids may be accepted if the following conditions are met:

- The device must be prepared to prevent accidental activation
- Battery terminals must be protected against short circuits
- The battery must either:
 - Remain securely attached and protected within the device, or
 - Be removed according to the manufacturer's instructions
- If the battery is removed:
 - It must be placed in **strong, rigid packaging**
 - It must be transported in the **cargo compartment only**
- Spare batteries:
 - No spare batteries allowed

⚠ Important Information

- Airline approval is required for all battery-powered wheelchairs and mobility aids
- Please contact the airline in advance of your travel to arrange approval and handling
- Always follow manufacturer instructions when preparing your device

