

Patient Sticker:



OIC Clavicle Plating System Inventory Control Sheet



Superior Distal Plates					
QTY	Description	Part #	GUDID #	Lot Number	Barcode
	4H LEFT (78mm)	CPD-04078L	10811998033829		(01)10811998033829
	6H LEFT (101mm)	CPD-06101L	10811998033843		(01)10811998033843
	8H LEFT (126mm)	CPD-08126L	10811998035403		(01)10811998035403
	10H LEFT (151mm)	CPD-10151L	10811998035427		(01)10811998035427
	4H RIGHT (78mm)	CPD-04078R	10811998033836		(01)10811998033836
	6H RIGHT (101mm)	CPD-06101R	10811998033850		(01)10811998033850
	8H RIGHT (126mm)	CPD-08126R	10811998035410		(01)10811998035410
	10H RIGHT (151mm)	CPD-10151R	10811998035434		(01)10811998035434



Midshaft Plates					
QTY	Description	Part #	GUDID #	Lot Number	Barcode
	8H LEFT (100mm)	CPM-08100L	10811998033874		(01)10811998033874
	10H LEFT (123mm)	CPM-10123L	10811998033898		(01)10811998033898
	8H RIGHT (100mm)	CPM-08100R	10811998033881		(01)10811998033881
	10H RIGHT (123mm)	CPM-10123R	10811998033904		(01)10811998033904
	6H NEUTRAL (76mm)	CPM-06076N	10811998033867		(01)10811998033867



Hole Cover					
QTY	Description	Part #	GUDID #	Lot Number	Barcode
	Hole Cover	LS-3500	10811998034468		(01)10811998034468

System Tips:

Trays Needed: OIC Small Fragment Base Tray, OIC Clavicle Plates

All 2.5mm Screws: Gold drill bits, Gold screw driver, Small depth gage

All 3.5mm Screws: White drill guides, White ring drill bits, White ring screwdriver, Large depth gage

Plates: Lime color for LEFT, Rose color for RIGHT, Gold color for NEUTRAL

OIC Customer Service: 800-619-2797

Fax: 775-636-8284

www.orthoimplantcompany.com