

Patient Sticker:



ORTEMIZ DRPx Distal Radius Plating
System Inventory Control Sheet



Narrow Distal Radius Plates					
QTY	Description	Part #	GUDID #	Lot Number	Barcode
	3H LEFT (45mm)	DRPX-N0345L	10811998035663		(01)10811998035663
	4H LEFT (53mm)	DRPX-N0453L	10811998035687		(01)10811998035687
	3H RIGHT (45mm)	DRPX-N0345R	10811998035670		(01)10811998035670
	4H RIGHT (53mm)	DRPX-N0453R	10811998035694		(01)10811998035694



Standard Distal Radius Plates					
QTY	Description	Part #	GUDID #	Lot Number	Barcode
	3H LEFT (51mm)	DRPX-S0351L	10811998035700		(01)10811998035700
	4H LEFT (59mm)	DRPX-S0459L	10811998035724		(01)10811998035724
	7H LEFT (89mm)	DRPX-S0789L	10811998035748		(01)10811998035748
	3H RIGHT (51mm)	DRPX-S0351R	10811998035717		(01)10811998035717
	4H RIGHT (59mm)	DRPX-S0459R	10811998035731		(01)10811998035731
	7H RIGHT (89mm)	DRPX-S0789R	10811998035755		(01)10811998035755



Dorsal Spanning Plates					
QTY	Description	Part #	GUDID #	Lot Number	Barcode
	8H (160mm)	DSP-08160	10811998039514		(01)10811998039514
	8H (210mm)	DSP-08210	10811998039521		(01)10811998039521

Disposables - PLEASE RE-ORDER

QTY	Description	Part #	GUDID #	Lot Number	Barcode
	1.6mm x 150mm	KW160-150	20811998034120		(01)20811998034120
	2.0mm x 95mm	DR20-095C	10811998039470		(01)10811998039470
	2.8mm x 95mm	DR28-095C	10811998039739		(01)10811998039739
	2.8mm x 155mm	DR28-155C	10811998034017		(01)10811998034017
	Short Guides	SG-0020	10811998039456		(01)10811998039456
	T6 Driver	01-T06090	10811998035632		(01)10811998035632

System Tips:

Trays Needed: OIC DRPx Distal Radius Plating System tray

Plates: Lime color for LEFT, Rose color for RIGHT

All Distal Screws: Purple ring drill guides, Purple ring drill bits, Purple ring screw driver


All Shaft Screws: White ring drill guides, White ring drill bits, White ring screwdriver



2.0 MM Locking Pegs (LPX-20XX)				
QTY	Length	Part #	GUDID #	Barcode
	14mm	LPX-2014	10811998035793	 (01)10811998035793
	16mm	LPX-2016	10811998035816	 (01)10811998035816
	18mm	LPX-2018	10811998035830	 (01)10811998035830
	20mm	LPX-2020	10811998035854	 (01)10811998035854
	22mm	LPX-2022	10811998035878	 (01)10811998035878
	24mm	LPX-2024	10811998035892	 (01)10811998035892
	26mm	LPX-2026	10811998035908	 (01)10811998035908
	28mm	LPX-2028	10811998035915	 (01)10811998035915



2.3 MM Locking Screws (LSX-23XX)				
QTY	Length	Part #	GUDID #	Barcode
	14mm	LSX-2314	10811998035953	 (01)10811998035953
	16mm	LSX-2316	10811998035977	 (01)10811998035977
	18mm	LSX-2318	10811998035991	 (01)10811998035991
	20mm	LSX-2320	10811998039012	 (01)10811998039012
	22mm	LSX-2322	10811998039036	 (01)10811998039036
	24mm	LSX-2324	10811998039050	 (01)10811998039050
	26mm	LSX-2326	10811998039067	 (01)10811998039067
	28mm	LSX-2328	10811998039074	 (01)10811998039074



2.6 MM Non-Locking Cortex Screws (STX-26XX)				
QTY	Length	Part #	GUDID #	Barcode
	8mm	STX-2608	10811998039241	 (01)10811998039241
	10mm	STX-2610	10811998039258	 (01)10811998039258
	12mm	STX-2612	10811998039265	 (01)10811998039265
	14mm	STX-2614	10811998039272	 (01)10811998039272
	16mm	STX-2616	10811998039296	 (01)10811998039296
	18mm	STX-2618	10811998039319	 (01)10811998039319
	20mm	STX-2620	10811998039333	 (01)10811998039333
	22mm	STX-2622	10811998039357	 (01)10811998039357
	24mm	STX-2624	10811998039371	 (01)10811998039371
	26mm	STX-2626	10811998039388	 (01)10811998039388
	28mm	STX-2628	10811998039395	 (01)10811998039395



2.6 MM Variable Angle Locking Screws (VLX-26XX)				
QTY	Length	Part #	GUDID #	Barcode
	8mm	VLX-2608	10811998039081	 (01)10811998039081
	10mm	VLX-2610	10811998039098	 (01)10811998039098
	12mm	VLX-2612	10811998039104	 (01)10811998039104
	14mm	VLX-2614	10811998039111	 (01)10811998039111
	16mm	VLX-2616	10811998039135	 (01)10811998039135
	18mm	VLX-2618	10811998039159	 (01)10811998039159
	20mm	VLX-2620	10811998039173	 (01)10811998039173
	22mm	VLX-2622	10811998039197	 (01)10811998039197
	24mm	VLX-2624	10811998039210	 (01)10811998039210
	26mm	VLX-2626	10811998039227	 (01)10811998039227
	28mm	VLX-2628	10811998039234	 (01)10811998039234



3.5 MM Non-Locking Cortex Screws (ST-35XX)				
QTY	Length	Part #	GUDID #	Barcode
	10mm	ST-3510	10811998035144	 (01)10811998035144
	11mm	ST-3511	10811998035588	 (01)10811998035588
	12mm	ST-3512	10811998035151	 (01)10811998035151
	13mm	ST-3513	10811998035595	 (01)10811998035595
	14mm	ST-3514	10811998035168	 (01)10811998035168
	15mm	ST-3515	10811998035601	 (01)10811998035601
	16mm	ST-3516	10811998035175	 (01)10811998035175
	18mm	ST-3518	10811998035182	 (01)10811998035182



3.5 MM Locking Screws (LS-35XX)				
QTY	Length	Part #	GUDID #	Barcode
	10mm	LS-3510	10811998034475	 (01)10811998034475
	12mm	LS-3512	10811998034482	 (01)10811998034482
	14mm	LS-3514	10811998034499	 (01)10811998034499
	16mm	LS-3516	10811998034505	 (01)10811998034505
	18mm	LS-3518	10811998034512	 (01)10811998034512