

Patient Sticker:



**ORTEMIZ Olecranon Plating
System Inventory Control Sheet**



Olecranon Plates					
QTY	Description	Part #	GUDID #	Lot Number	Barcode
	3H LEFT (78mm)	OLEC-03078L	10811998039784		 (01)10811998039784
	6H LEFT (108mm)	OLEC-06108L	10811998039807		 (01)10811998039807
	10H LEFT (152mm)	OLEC-10152L	10811998039821		 (01)10811998039821
	3H RIGHT (78mm)	OLEC-03078R	10811998039791		 (01)10811998039791
	6H RIGHT (108mm)	OLEC-06108R	10811998039814		 (01)10811998039814
	10H RIGHT (152mm)	OLEC-10152R	10811998039838		 (01)10811998039838

System Tips:

Trays Needed: OIC Small Fragment Base Tray, OIC Olecranon Plates

All 3.5mm Screws: White drill guides, White ring drill bits, White ring screwdriver, Large depth gage

Plates: Lime color for LEFT, Rose color for RIGHT