

Patient Sticker:



ORTEMIZ Proximal Humerus Plating  
System Inventory Control Sheet



| Proximal Humerus Plates |                  |            |                |            |                                                                                                        |
|-------------------------|------------------|------------|----------------|------------|--------------------------------------------------------------------------------------------------------|
| QTY                     | Description      | Part #     | GUDID #        | Lot Number | Barcode                                                                                                |
|                         | 3H LEFT (100mm)  | PHP-03100L | 10811998034840 |            |  (01)10811998034840 |
|                         | 5H LEFT (124mm)  | PHP-05124L | 10811998034864 |            |  (01)10811998034864 |
|                         | 7H LEFT (148MM)  | PHP-07148L | 10811998034888 |            |  (01)10811998034888 |
|                         | 9H LEFT (172MM)  | PHP-09172L | 10811998039845 |            |  (01)10811998039845 |
|                         | 3H RIGHT (100mm) | PHP-03100R | 10811998034857 |            |  (01)10811998034857 |
|                         | 5H RIGHT (124mm) | PHP-05124R | 10811998034871 |            |  (01)10811998034871 |
|                         | 7H RIGHT (148MM) | PHP-07148R | 10811998034895 |            |  (01)10811998034895 |
|                         | 9H RIGHT (172MM) | PHP-09172R | 10811998039852 |            |  (01)10811998039852 |

System Tips:

Trays Needed: OIC Small Fragment Base Tray, OIC Proximal Humerus Plates

All 3.5mm Screws: White drill guides, White ring drill bits, White ring screwdriver, Large depth gage

Plates: Lime color for LEFT, Rose color for RIGHT