

Patient Sticker:



ORTEMIZ Clavicle X Plating

System Inventory Control Sheet



Lateral Anterior Clavicle Plates					
QTY	Description	Part #	GUDID #	Lot Number	Barcode
	8H NEUTRAL (69mm)	CAL-08069	10811998031948		 (01)10811998031948
	9H NEUTRAL (81mm)	CAL-09081	10811998031955		 (01)10811998031955



Midshaft Anterior Plates					
QTY	Description	Part #	GUDID #	Lot Number	Barcode
	6H NEUTRAL (71mm)	CAM-06071	10811998031962		 (01)10811998031962
	8H NEUTRAL (92mm)	CAM-08092	10811998031979		 (01)10811998031979
	10H NEUTRAL (113mm)	CAM-10113	10811998031986		 (01)10811998031986



Hole Cover					
QTY	Description	Part #	GUDID #	Lot Number	Barcode
	Hole Cover	LS-3500	10811998034468		 (01)10811998034468

System Tips:

Trays Needed: OIC Small Fragment Base Tray, OIC Clavicle X Plates

All 2.5mm Screws: Gold drill bits, Gold screw driver, Small depth gage

All 3.5mm Screws: White drill guides, White ring drill bits, White ring screwdriver, Large depth gage

Plates: All plates can be used on either right or left anatomy