

Patient Sticker:



ORTEMIZ Ankle Plating System
Inventory Control Sheet



Distal Fibula Plates					
QTY	Description	Part #	GUDID #	Lot Number	Barcode
	3H NEUTRAL (65mm)	DFIB-03065	10811998033911		(01)10811998033911
	4H NEUTRAL (77mm)	DFIB-04077	10811998033928		(01)10811998033928
	5H NEUTRAL (87mm)	DFIB-05087	10811998033935		(01)10811998033935
	7H NEUTRAL (109mm)	DFIB-07109	10811998033942		(01)10811998033942
	9H NEUTRAL (131mm)	DFIB-09131	10811998033959		(01)10811998033959
	11H NEUTRAL (153mm)	DFIB-11153	10811998033966		(01)10811998033966



Medial Distal Tibia Plates					
QTY	Description	Part #	GUDID #	Lot Number	Barcode
	3H RIGHT (89mm)	MDT-03089R	10811998034734		(01)10811998034734
	6H RIGHT (127mm)	MDT-06127R	10811998034758		(01)10811998034758
	8H RIGHT (152mm)	MDT-08152R	10811998034772		(01)10811998034772
	10H RIGHT (177mm)	MDT-10177R	10811998034796		(01)10811998034796
	13H RIGHT (214mm)	MDT-13214R	10811998034819		(01)10811998034819
	16H RIGHT (252mm)	MDT-16252R	10811998034833		(01)10811998034833
	3H LEFT (89mm)	MDT-03089L	10811998034727		(01)10811998034727
	6H LEFT (127mm)	MDT-06127L	10811998034741		(01)10811998034741
	8H LEFT (152mm)	MDT-08152L	10811998034765		(01)10811998034765
	10H LEFT (177mm)	MDT-10177L	10811998034789		(01)10811998034789
	13H LEFT (214mm)	MDT-13214L	10811998034802		(01)10811998034802
	16H LEFT (252mm)	MDT-16252L	10811998034826		(01)10811998034826



Anterolateral Distal Tibia Plates					
QTY	Description	Part #	GUDID #	Lot Number	Barcode
	4H RIGHT (69mm)	ALDT-04069R	10811998033652		 (01)10811998033652
	6H RIGHT (95mm)	ALDT-06095R	10811998033676		 (01)10811998033676
	8H RIGHT (120mm)	ALDT-08120R	10811998033690		 (01)10811998033690
	10H RIGHT (145mm)	ALDT-10145R	10811998033713		 (01)10811998033713
	13H RIGHT (128mm)	ALDT-13128R	10811998033737		 (01)10811998033737
	4H LEFT (69mm)	ALDT-04069L	10811998033645		 (01)10811998033645
	6H LEFT (95mm)	ALDT-06095L	10811998033669		 (01)10811998033669
	8H LEFT (120mm)	ALDT-08120L	10811998033683		 (01)10811998033683
	10H LEFT (145mm)	ALDT-10145L	10811998033706		 (01)10811998033706
	13H LEFT (128mm)	ALDT-13128L	10811998033720		 (01)10811998033720

System Tips:

Trays Needed: ORTEMIZ Small Fragment Base Tray, OIC Ankle Plates

All 2.5mm Screws: Gold drill bits, Gold screw driver, Small depth gage

All 3.5mm Screws: White drill guides, White ring drill bits, White ring screwdriver, Large depth gage

Plates: Lime color for LEFT, Rose color for RIGHT, Gold color for NEUTRAL