

Patient Sticker:



ORTEMIZ ARFN Nail

System Inventory

Control Sheet



| Antegrade/Retrograde Femur Nails |                |             |                |            |                                                                                                          |
|----------------------------------|----------------|-------------|----------------|------------|----------------------------------------------------------------------------------------------------------|
| QTY                              | Description    | Part #      | GUDID #        | Lot Number | Barcode                                                                                                  |
|                                  | 10mm x 320mm   | ARFN-100320 | 10811998036608 |            |  (01)10811998036608   |
|                                  | 11.5mm x 320mm | ARFN-115320 | 10811998036677 |            |  (01)10811998036677   |
|                                  | 13mm x 320mm   | ARFN-130320 | 10811998036745 |            |  (01)10811998036745   |
|                                  | 10mm x 340mm   | ARFN-100340 | 10811998036615 |            |  (01)10811998036615   |
|                                  | 11.5mm x 340mm | ARFN-115340 | 10811998036684 |            |  (01)10811998036684   |
|                                  | 13mm x 340mm   | ARFN-130340 | 10811998036752 |            |  (01)10811998036752   |
|                                  | 10mm x 360mm   | ARFN-100360 | 10811998036622 |            |  (01)10811998036622   |
|                                  | 11.5mm x 360mm | ARFN-115360 | 10811998036691 |            |  (01)10811998036691   |
|                                  | 13mm x 360mm   | ARFN-130360 | 10811998036769 |            |  (01)10811998036769   |
|                                  | 10mm x 380mm   | ARFN-100380 | 10811998036639 |            |  (01)10811998036639 |
|                                  | 11.5mm x 380mm | ARFN-115380 | 10811998036707 |            |  (01)10811998036707 |
|                                  | 13mm x 380mm   | ARFN-130380 | 10811998036776 |            |  (01)10811998036776 |
|                                  | 10mm x 400mm   | ARFN-100400 | 10811998036646 |            |  (01)10811998036646 |
|                                  | 11.5mm x 400mm | ARFN-115400 | 10811998036714 |            |  (01)10811998036714 |
|                                  | 13mm x 400mm   | ARFN-130400 | 10811998036783 |            |  (01)10811998036783 |
|                                  | 10mm x 420mm   | ARFN-100420 | 10811998036653 |            |  (01)10811998036653 |
|                                  | 11.5mm x 420mm | ARFN-115420 | 10811998036721 |            |  (01)10811998036721 |
|                                  | 13mm x 420mm   | ARFN-130420 | 10811998036790 |            |  (01)10811998036790 |
|                                  | 10mm x 440mm   | ARFN-100440 | 10811998036660 |            |  (01)10811998036660 |
|                                  | 11.5mm x 440mm | ARFN-115440 | 10811998036738 |            |  (01)10811998036738 |
|                                  | 13mm x 440mm   | ARFN-130440 | 10811998036806 |            |  (01)10811998036806 |



| 5.0mm Cross-Lock Screws |        |           |                |            |                                                                                                          |
|-------------------------|--------|-----------|----------------|------------|----------------------------------------------------------------------------------------------------------|
| QTY                     | Length | Part #    | GUDID #        | Lot Number | Barcode                                                                                                  |
|                         | 25mm   | NLB-5025  | 10811998037117 |            |  (01)10811998037117   |
|                         | 27.5mm | NLB-5027  | 10811998037124 |            |  (01)10811998037124   |
|                         | 30mm   | NLB-5030  | 10811998037131 |            |  (01)10811998037131   |
|                         | 32.5mm | NLB-5032  | 10811998037148 |            |  (01)10811998037148   |
|                         | 35mm   | NLB-5035  | 10811998037155 |            |  (01)10811998037155   |
|                         | 37.5mm | NLB-5037  | 10811998037162 |            |  (01)10811998037162   |
|                         | 40mm   | NLB-5040  | 10811998037179 |            |  (01)10811998037179   |
|                         | 42.5mm | NLB-5042  | 10811998037186 |            |  (01)10811998037186   |
|                         | 45mm   | NLB-5045  | 10811998037193 |            |  (01)10811998037193   |
|                         | 47.5mm | NLB-5047  | 10811998037209 |            |  (01)10811998037209   |
|                         | 50mm   | NLB-5050  | 10811998037216 |            |  (01)10811998037216   |
|                         | 55mm   | NLB-5055  | 10811998037223 |            |  (01)10811998037223   |
|                         | 60mm   | NLB-5060  | 10811998037230 |            |  (01)10811998037230   |
|                         | 65mm   | NLB-5065  | 10811998037247 |            |  (01)10811998037247   |
|                         | 70mm   | NLB-5070  | 10811998037254 |            |  (01)10811998037254   |
|                         | 75mm   | NLB-5075  | 10811998037261 |            |  (01)10811998037261   |
|                         | 80mm   | NLB-5080  | 10811998037278 |            |  (01)10811998037278  |
|                         | 85mm   | NLB-5085  | 10811998037285 |            |  (01)10811998037285 |
|                         | 90mm   | NLB-5090  | 10811998037292 |            |  (01)10811998037292 |
|                         | 95mm   | NLB-5095  | 10811998037308 |            |  (01)10811998037308 |
|                         | 100mm  | NLB-50100 | 10811998037315 |            |  (01)10811998037315 |

| Disposables - PLEASE RE-ORDER |                         |             |                |            |                                                                                                          |
|-------------------------------|-------------------------|-------------|----------------|------------|----------------------------------------------------------------------------------------------------------|
| QTY                           | Description             | Part #      | GUDID #        | Lot Number | Barcode                                                                                                  |
|                               | 3.2mm Guidepin (Each)   | 32-400      | 10811998036592 |            |  (01)10811998036592 |
|                               | 3.2mm Guidepin (5-Pack) | 32-400P     | 20811998036599 |            |  (01)20811998036599 |
|                               | Ball Tipped Guiderod    | GW03-1000-S | 10811998036851 |            |  (01)10811998036851 |
|                               | 4.2mm x 330mm Drillbit  | DR42-330C   | 10811998036837 |            |  (01)10811998036837 |
|                               | 4.2mm x 155mm Drillbit  | DR42-155    | 10811998036813 |            |  (01)10811998036813 |

#### System Tips:

Trays Needed: OIC IM Nail Base Tray, OIC Femur/Hip Tray, ARFN Implant Tray

Use 12.5mm Reamer for Entry

Nail can be used Retrograde or Antegrade (periformis)