

| MRI IAC W & WO       |                  |          |                      |                 |
|----------------------|------------------|----------|----------------------|-----------------|
| Technical Parameters |                  |          |                      |                 |
| Anatomic Coverage:   | Centered on IACs |          |                      |                 |
| Coil:                | Head             |          |                      |                 |
| Scanner Preference:  |                  |          |                      |                 |
| Contrast:            | Pre-injection:   |          |                      | Post-injection: |
| Scanning Plane:      | Axial            | Axial    | Axial                | Axial           |
| Slice Thickness:     | 0.5mm            | 0.5mm    | 2mm                  | 0.5mm           |
| Gap:                 | 0mm              | 0mm      | 0.2mm                | 0mm             |
| Weighting:           | T2               | T1       | DWI                  | T1              |
| Mode:                | 3D TSE           | 3D TSE   | 2D TSE or<br>RESOLVE | 3D TSE          |
| FAT SAT:             | NO               | NO       | YES                  | YES             |
| FOV:                 | 17cm             | 17cm     | 17cm                 | 17cm            |
| MPR:                 | COR, SAG, OBLIQ  | COR, SAG | NO                   | COR, SAG        |
| Notes:               |                  |          |                      |                 |