

MRI Neck Soft Tissue W & WO							
Technical Parameters							
Anatomic Coverage:	Mid orbits to mid aortic arch.						
Coil:	NV or Spine						
Scanner Preference:							
Contrast:	Pre-Injection:				Post-injection:		
Scanning Plane:	Sagittal	Axial	Coronal	Axial	Sagittal	Coronal	Axial
Slice Thickness:	0.7mm	3mm	2.5mm	3mm	0.7mm	2mm	2mm
Gap:	0mm	0.3mm	0.25mm	0.3mm	0mm	0.25mm	0mm
Weighting:	T1	DWI	T2	T2	T1	T1	T1
Mode:	3D TSE	2D EPI	2D TSE	2D TSE	3D TSE	2D TSE	3D VIBE/THRIVE
FAT SAT:	NO	YES	DIXON: IN PHASE	DIXON: IN PHASE	YES	DIXON: IN PHASE	YES
FOV:	21cm	22cm	22cm	21cm	21cm	22cm	18cm
MPR:	NO	NO	NO	NO	AXIAL, COR	NO	NO
Notes:							