

MRI TMJ					
Technical Parameters					
Anatomic Coverage:	Skull base to vertex				
Coil:	Head or NV				
Scanner Preference:					
Scanning Plane:	Axial	Sagittal (R to L)	Sagittal (R to L)	Coronal (R to L)	Sagittal (R to L)
Slice Thickness:	4mm	2mm	2mm	2mm	2mm
Gap:	1mm	0mm	0mm	0mm	0mm
Weighting:	T1	PD Closed	T2 Closed	T1 Closed	PD Open
Mode:	TSE	TSE	TSE	TSE	TSE
FAT SAT:	NO	NO	NO	NO	NO
FOV:	23cm	12cm	12cm	12cm	12cm
MPR:	NO	NO	NO	NO	NO
Notes:					