

MRI Bony Pelvis Fast Fracture WO		
Technical Parameters		
Anatomic Coverage:	Entire pelvis	
Coil:	Torso	
Scanner Preference:	1.5T or 3T	
Scanning Plane:	Coronal	Coronal
Slice Thickness:	5-6mm	5-6mm
Gap:	1mm	1mm
Weighting:	T1	T2 SPIR
Mode:	TSE	TSE
FAT SAT:	NO	STIR
FOV:	30cm	30cm
MPR:	NO	NO
Notes:		