

| MRI Forearm W & WO   |                            |                 |         |          |
|----------------------|----------------------------|-----------------|---------|----------|
| Technical Parameters |                            |                 |         |          |
| Anatomic Coverage:   | Elbow through wrist joints |                 |         |          |
| Coil:                | XL Torso                   |                 |         |          |
| Scanner Preference:  | 1.5T                       |                 |         |          |
| Contrast:            | Pre-injection:             | Post-injection: |         |          |
| Scanning Plane:      | Axial                      | Axial           | Coronal | Sagittal |
| Slice Thickness:     | 4mm                        | 4mm             | 4mm     | 4mm      |
| Gap:                 | 1mm                        | 1mm             | 1mm     | 1mm      |
| Weighting:           | T1 FS                      | T1 FS           | T1 FS   | T1 FS    |
| Mode:                | TSE                        | TSE             | TSE     | TSE      |
| FAT SAT:             | YES                        | YES             | YES     | YES      |
| FOV:                 | 14cm                       | 14cm            | 14cm    | 14cm     |
| MPR:                 | NO                         | NO              | NO      | NO       |
| Notes:               |                            |                 |         |          |