

MRI Osteomyelitis						
Technical Parameters						
Anatomic Coverage:	Area of interest					
Coil:	As appropriate for anatomic location					
Scanner Preference:	1.5T or 3T					
Scanning Plane:	Axial	Axial	Sagittal	Sagittal	Coronal	Coronal
Slice Thickness:	As appropriate for anatomic location					
Gap:						
Weighting:	T1	STIR	T1	STIR	T1	STIR
Mode:	TSE	TSE	TSE	TSE	TSE	TSE
FAT SAT:	NO	YES	NO	YES	NO	YES
FOV:	As appropriate for anatomic location					
MPR:	NO	NO	NO	NO	NO	NO
Notes:						