

MRI Sacrum/SI Joints W & WO									
Technical Parameters									
Anatomic Coverage:	SI joints, sacrum, coccyx								
Coil:	Torso								
Scanner Preference:	1.5T or 3T								
Contrast:	Pre-injection:						Post-injection:		
Scanning Plane:	Axial T1	Axial PDFS	Sagittal T2FS	Coronal Oblique T1	Coronal Oblique PDFS	Coronal STIR	Coronal Oblique	Axial Oblique	Sagittal Oblique
Slice Thickness:	3 mm	3 mm	3 mm	3 mm	3 mm	4 mm	3mm	3mm	3mm
Gap:	0.5 mm	0.5 mm	0.5 mm	0.5 mm	0.5mm	1 mm	0.5mm	0.5mm	0.5mm
Weighting:	T1	PD Fat-Sat	T2 Fat-Sat	T1	PD Fat-Sat	STIR	T1 FS	T1 FS	T1 FS
Mode:	TSE	TSE	TSE	TSE	TSE	TSE	TSE	TSE	TSE
FAT SAT:	NO	YES	YES	NO	YES	YES	YES	YES	YES
FOV:	18-22 cm	18-22 cm	18-22 cm	18-22 cm	18-22 cm	30-40 cm	18-22 cm	18-22 cm	18-22 cm
MPR:	NO	NO	NO	NO	NO	NO	NO	NO	
Notes:									

