

MRI Soft Tissue Tumor

Technical Parameters

Anatomic Coverage:	Area of interest							
Coil:	As appropriate for anatomic location							
Scanner Preference:	1.5T or 3T							
Contrast:	Pre-contrast:					Post-contrast:		
Scanning Plane:	Axial	Axial	Sagittal	Coronal	Axial	Axial	Sagittal	Coronal
Slice Thickness:	As appropriate for anatomic location					As appropriate for anatomic location		
Gap:								
Weighting:	T1	T2	T2	T2	T1	T1	T1	T1
Mode:	TSE	TSE	TSE	TSE	TSE	TSE	TSE	TSE
FAT SAT:	NO	YES	NO	YES	YES	YES	YES	YES
FOV:	As appropriate for anatomic location					As appropriate for anatomic location		
MPR:	NO	NO	NO	NO	NO	NO	NO	NO
Notes:								