

| MRI Tib/Fib W & WO   |                 |                 |         |          |
|----------------------|-----------------|-----------------|---------|----------|
| Technical Parameters |                 |                 |         |          |
| Anatomic Coverage:   | Knees to Ankles |                 |         |          |
| Coil:                | Torso           |                 |         |          |
| Scanner Preference:  | 1.5T or 3T      |                 |         |          |
| Contrast:            | Pre-injection:  | Post-injection: |         |          |
| Scanning Plane:      | Axial           | Axial           | Coronal | Sagittal |
| Slice Thickness:     | 5mm             | 5mm             | 5mm     | 5mm      |
| Gap:                 | 1mm             | 1mm             | 1mm     | 1mm      |
| Weighting:           | T1 FS           | T1 FS           | T1 FS   | T1 FS    |
| Mode:                | TSE             | TSE             | TSE     | TSE      |
| FAT SAT:             | YES             | YES             | YES     | YES      |
| FOV:                 | 40cm            | 40cm            | 40cm    | 40cm     |
| MPR:                 | NO              | NO              | NO      | NO       |
| Notes:               |                 |                 |         |          |